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# NEBULIZER Aerosol Therapy System

## INSTRUCTION MANUAL

Thank you for choosing our Nebulizer Gamma® designed to carry out aerosol therapy of acute and chronic diseases of upper and lower respiratory passages. This device allows for easy and effective transportation of big medicine doses directly to the disease centre.

We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Gamma Trade Mark.

Before starting to use the nebulizer system please study the instruction manual carefully. The instruction manual offers all information you need to correctly carry out aerosol therapy.

### ATTENTION

Please use only those medicines that were prescribed to you by your doctor, DO NOT execute self-treatment and DO NOT change the methods prescribed without consulting your doctor.


## CONTENT

1. COMPONENTS AND DESCRIPTION OF YOUR NEBULIZER SYSTEM  
MODEL CHU-CHU-TRAIN
2. GETTING READY FOR USAGE
  - 2.1. ASSEMBLING THE DEVICE
  - 2.2. REPLACEMENT OF FILTER
3. CARRYING OUT AEROSOL THERAPY
4. CLEANING AND DISINFECTION
5. POSSIBLE PROBLEMS AND HOW TO SOLVE THEM
6. TECHNICAL SPECIFICATIONS
7. WARRANTEE
8. WHAT IS NEBULIZING THERAPY

## 1. COMPONENTS AND DESCRIPTION OF YOUR NEBULIZER MODEL CHU-CHU-TRAIN

Your nebulizer system includes compressor unit Chu-Chu-Train, power cord, PVC tubing, mouth-piece, jet nebulizer, 3 air filters, adult mask, pediatric mask, Instruction Book, guarantee card, packing.

The picture shows nebulizer model Chu-Chu-Train

### ATTENTION:

Read the manual carefully before using the device and keep it in case you need further information in the future.

This nebulizer must be used only for the specific cause it was created for: the aerosoltherapy. Any other use will be considered unsuitable and thus dangerous; the manufacturer can not be considered responsible for the damages caused by unsuitable, incorrect or irrational use, or if the unit has been used in electrical installations not in accordance with the standing safety standards.

- When using any electrical appliance certain important safety measures must always be observed, including the following:
  - The unit is not protected from liquid penetration.
  - The nebulizer system must not be touched with wet or soaked hands.
  - Do not use the device in presence of inflammable anesthetic mix in contact with the air, oxygen or nitrogen protoxide.
  - Do not block the air ventilation holes.
  - Unplug the unit before cleaning or maintaining it.
  - Disconnect the nebulizer system when it is not in use.
- Do not untwist the jet nebulizer when the unit is operating.
- Do not leave medication inside the jet nebulizer when it is not in use.
- Do not pour more than 10 ml in the jet nebulizer.
- Use the unit on stable surfaces. Never use it while sleepy or drowsy.
- Use the unit intermittently. Maximum operating period: 30 minutes. It must stand unused for 30 minutes.
- The mask and the mouth piece are for personal use, and they have not to be shared in order to avoid any infection or contagion.
- Stop the nebulization if you feel any indisposition.
- Children will use this unit only under adult supervision.
- This unit have to be repaired only by authorized service centres.

**NOTE:** this unit employs a thermal limiter that will shut off the entire product should the motor overheat.

## 2. GETTING READY FOR USAGE

### 2.1. ASSEMBLING THE DEVICE

1. Take the unit out of the box.

2. Open the jet nebulizer (Fig. A) and pour the medication to nebulize into it.

Do not pour more than 10 ml (Fig. B).


Then twist both parts and close again. (Fig. C).

3. Connect the tubing to the jet nebulizer and the other end to the nebulizer outlet. (Fig. D).

4. Put the mask or the mouth piece to be used into the jet nebulizer (Fig. E).

5. Connect the unit to a suitable power receptacle.

Your nebulizer system is ready for work.

### 2.2. REPLACEMENT OF THE FILTER

1. Filter should be changed every month or sooner if the filter turns completely gray in color.

2. Remove filter cap by grasping it firmly and pulling out the unit.

3. Change the dirty filter of the cap with small, pointed object.

4. Insert a clean San Up filter. Extra filters are available.

5. Push the filter cap with new filter into position.

### Caution:

Do not use cotton or any other material. Do not wash or clean the filter. Use only original filters. Do not operate without a filter.

## 3. CARRYING OUT AEROSOL THERAPY

1. Turn the unit by pressing the ON/OFF switch to start the treatment.

### ATTENTION:

This device is for intermittent use: 30 min. ON/30 min. OFF.

2. Inhale the aerosol solution using the prescribed accessory in calm way. One should sit at the table with one's back upright, relaxed; otherwise you will constrict your respiratory passages and impair the effect of inhalation.

3. Turn the unit off and unplug it once the treatment has finished.

4. Clean the jet nebulizer and the devices as indicated in Cleaning and disinfection section. Store the jet nebulizer and the devices into the unit's lodging.

## 4. CLEANING AND DISINFECTION

### Important:

1. Wash all the pieces, except for the air filter, with warm water and detergent. Rinse with warm tap water to remove the detergent.

2. Dip the pieces into a clean bowl containing 3 parts of warm water and 1 part of vinegar during 30 minutes, disinfectant can also be used.

3. Wash your hands and take out the pieces from the disinfectant solution, rinse them with warm tap water and leave them on an absorbent piece of paper.

4. Keep the pieces into a hermetic bag. Do not dry them with a towel.

5. Clean the outer surface with a dry piece of cloth.

- Do not use any volatile fluid like alcohol or benzene.

- The jet nebulizer and the devices must NOT be boiled.

- The jet nebulizer has to be replaced if it is stuck by medication or other substances.

- Under no circumstances will the unit be opened.

## 5. POSSIBLE PROBLEMS AND HOW TO SOLVE THEM

If something wrong happens while the device is working, please check and correct the following points:

Symptom	Possible cause	Suggested solutions
Although the device is on, it does not work	1. Incorrectly plugged 2. No power	1. Check the plug is properly connected 2. Check the power is not gone out
The mist does not come out	1. There is no solution in the jet nebulizer 2. There is too much solution in the jet nebulizer 3. PVC tubing is disconnected	1. Add medication, maximum 10 ml 2. Remove excess medication from cup. 3. Connect the PVC tubing.
Aerosol output reduced (low)	1. Air filter is clogged. 2. The jet nebulizer is dirty.	1. Replace for a clean air filter. 2. Replace for new one.

Should the device not start working properly again, please contact the Service Centre; the address is stated in the Guarantee Card.

## 6. TECHNICAL SPECIFICATIONS

Method of nebulization	Dry piston compressor
Operating Voltage	220V, 50Hz
Class	
Power consumption	0.6A
Jet nebulizer maximum capacity	10 ml

Sound noise level (at 1m distance)	60 dB
Mean Mass Aerodynamic Diameter (microns)	less than 5 microns
Medication outlet*	0.33 ml/min. or greater; value measured with 0.9% NaCl solution.
Weight	1.732 kg
Unit dimensions	240 x 140 x 124 mm
Operation Humidity and Temperature	-5 to +40° C, 30% to 95%
Storage Humidity and Temperature	-10 to +40° C, 10% to 95%
Atmospheric pressure range for carriage and storage	50 - 106 kPa (500-1060 mbar)

\*This specification may vary with the kind of the medication being used and its viscosity. The specifications and the outer design are subject to change for improvement reasons without previous notice.

### Symbols Key

	ON
	OFF
	Type B
	Read the instruction
	AC
	Class II

## 7. WARRANTY

The nebulizer system Chu-Chu-Train is covered by a 2-year warranty from date of purchase for any defects in the design or in the materials used. The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear.

The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.

The manufacturer shall not be liable for any damages to people, property or pets, either directly or indirectly, caused by the non-compliance with the instruction book.

The warrantee is only valid upon presentation of the warrantee card which was correctly filled in and sealed.

## 8. WHAT IS NEBULIZING THERAPY

### What is a nebulizer?

The name of the device comes from the Latin word nebula that means fog, cloud. The nebulizer system is a device that turns a liquid drug into a fine aerosol under the influence of the compressed air. Aerosol is a drug suspension. So, nebulizing therapy is an inhalation therapy with the help of aerosols.

Inhalation therapy has been used for treatment of lung diseases for a long time. The ancient civilizations of Egypt, India, China, the Middle East, Greece are known to have used inhalations of evaporations from menthol, eucalyptus, different plants.

For patients with lung diseases nebulizing therapy seems to be the most rational one since big drug doses come very fast to the respiratory tract directly.

### How does the nebulizer system work?

Nebulizer directs an aerosol containing larger and smaller particles to the respiratory tract. The gravity makes the larger particles accumulate on the walls of the upper respiratory passages while the smaller particles have enough time to reach the lower bronchial branches. Therefore the drug

comes to all parts of the patient's bronchial tree, including the smallest bronchi and alveole, and gets into the blood vessels on the walls of the alveole in big therapeutic doses with tiny or none side effects. Some drugs are very slow to be absorbed from the alimentary canal or are affected dramatically by the first passage through the liver. In such cases the way out is to use inhalation therapy. Nebulizing therapy can be used by patients of any age, it's suitable for small children and infirm.

### What is nebulizing therapy used for?

Nebulizer systems are used by the patients with:

- acute respiratory diseases;
- obstructive pulmonary disease;
- chronic inflammatory processes in respiratory passages such as chronic bronchitis, bronchial asthma (incl. chronic or unstable ones, exacerbation of asthma), chronic pharyngitis;
- respiratory passages infections;
- abnormalities in secretion quality and secreting in bronchi;
- different types of coughing;

Nebulizing therapy helps to decrease the level of inflammation activity, to improve respiratory passage drain function, to lessen the bronchismus, to make the edema smaller and to protect the mucous membrane, etc.

### What are the advantages of nebulizing therapy compared to other means of drug delivery?

Nebulizing therapy gives a possibility to quickly effect the inflammatory zone by big drug doses which increases the influence of the treatment and reduces the risk of side effects to other organs. Nebulizing therapy is the easiest way of inhalation treatment and it is also an economical means of treatments from the point of view of purposeful use of drug. Aerosol treatment with the help of nebulizing systems can be used for children, elderly and infirm patients.

**Please use only those drugs that were prescribed to you by your doctor, DO NOT execute self-treatment and DO NOT change the methods/frequency prescribed without consulting your doctor.**

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